




# Rosebud Reminder Clock User Guide



## To Get To The Settings Menu

Press 'MENU' button and it will display the MENU settings. Press 'MENU' again to exit. If no action is taken within 30 seconds, it will return to the main display automatically. Please note: Changes will be made to the settings highlighted .

## Set Time /Date

Select 'set time' or 'set date' as needed so that it is highlighted and press 'OK' . While the digit is underlined , press 'up' and 'down' to change the time. When you get to the correct digit press 'Left' and 'right' to move onto the next digit. Press 'OK' to confirm or time/date will not be changed.

## Set 12hr mode/ 24h mode

In settings menu use 'up' and 'down' buttons to scroll to 'time mode' . When this is highlighted use 'left' and 'right' to identify your required time mode either 12hr or 24hr, press 'OK' to confirm.

## Set Calendar Mode

To choose from our 4 display settings (calendar modes) go to the settings menu and scroll down to 'Calendar mode' When this is highlighted. Use 'left' and 'right' to select the calendar mode you require.

## Set Brightness - Day/ Night

In settings menu highlight 'brightness- day' or 'brightness- night'. Use 'left' and 'right' to make the adjustment. The larger digit will have a brighter luminance

## Set contrast

In settings menu highlight 'set contrast'. Use ' left' and 'right' buttons to make adjustments.

**Set one of Rosebuds pre programmed reminders (Up to 20 reminders can be set daily)**

in the settings Menu, scroll using 'up' and 'down' buttons to one of the alarms (01 - 20). When your chosen alarm is highlighted in red press 'OK' twice to bring up the list of reminders. Use 'up' and 'down' to highlight your required reminder and press 'OK' again. Press 'left' and right to underline the time and press 'up' and 'down' to select the time at which you wish the reminder to sound. Once the correct time is selected use 'left' and 'right' followed by 'up' and 'down' to select 'Repeat', 'once' or 'off'. ('Repeat' will play the reminder at your selected time everyday, 'Once' will play the reminder once only, 'Off' will not play the reminder). When you have selected the reminder and the time you require press 'OK' to confirm.

**Personalising your rosebud reminders****To add text and sounds reminders.**

First you must create your text or picture. The easiest way to do this is by using a program such as Microsoft Paint (you can use many other similar programs) Open up 'Paint' -File - Properties select image size 1024pixels wide X 768pixels High then create your chosen image and save. Create a folder in your chosen USB or SD card and drop the image into the folder. To add sound make a sound recording and drop into the same folder. You can create as many folders like this as your USB or SD card memory will hold. When this USB/ SD Card is plugged into your Rosebud Reminder Clock you can then select from your personalised list of reminders by the method set out above.

**To add video reminders**

Make sure your video recording is in landscape. Create a folder in your chosen USB or SD Card and drop the video into this folder. When this USB/SD card is then plugged into your Rosebud Reminder clock you can select from your personalised list of reminders as set out above.

The Rosebud Reminder Clock was designed in association with Health Care Professionals specifically for cognitive impairment, Dementia, Alzheimer's, Autism, Memory Loss, Brain Injury and Elderly care.

Rosebuds has a friendly voice and can be a helpful reminder and company for your loved one when you can't be there, reducing anxiety and stress and helping people to keep their independence for longer.

**Ravencourt Ltd**

**Drift Road**

**Stamford**

**PE9 1UT**

**[www.ravencourt.com](http://www.ravencourt.com)**

**[info@ravencourt.com](mailto:info@ravencourt.com)**

**01780 489100**